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
Specimen type	Collection materials	Transport to laboratory	Storage till testing	Comment
Nasopharyngeal and oropharyngeal swab	Dacron or polyester flocked swabs*	4 °C	≤5 days: 4 °C >5 days: -70 °C	The nasopharyngeal and oropharyngeal swabs should be placed in the same tube to increase the viral load.
Bronchoalveolar lavage	sterile container *	4 °C	≤48 hours: 4 °C >48 hours: -70 °C	There may be some dilution of pathogen, but still a worthwhile specimen
(Endo)tracheal aspirate, nasopharyngeal aspirate or nasal wash	sterile container *	4 °C	≤48 hours: 4 °C >48 hours: -70 °C	
Sputum	sterile container	4 °C	≤48 hours: 4 °C >48 hours: -70 °C	Ensure the material is from the lower respiratory tract
Tissue from biopsy or autopsy including from lung	sterile container with saline	4 °C	≤24 hours: 4 °C >24 hours: -70 °C	
Serum (2 samples acute and convalescent possibly 2-4 weeks after acute phase)	Serum separator tubes (adults: collect 3-5 ml whole blood)	4 °C	≤5 days: 4 °C >5 days: -70 °C	Collect paired samples: • acute – first week of illness • convalescent – 2 to 3 weeks later
Whole blood	collection tube	4 °C	≤5 days: 4 °C >5 days: -70 °C	For antigen detection particularly in the first week of illness
Urine	urine collection container	4 °C	≤5 days: 4 °C >5 days: -70 °C	

*For transport of samples for viral detection, use VTM (viral transport medium) containing antifungal and antibiotic supplements. For bacterial or fungal culture, transport dry or in a very small amount of sterile water. Avoid repeated freezing and thawing of specimens.

Aside from specific collection materials indicated in the table also assure other materials and equipment are available: e.g. transport containers and specimen collection bags and packaging, coolers and cold packs or dry ice, sterile blood-drawing equipment (e.g. needles, syringes and tubes), labels and permanent markers, PPE, materials for decontamination of surfaces.

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۱. C Moore, S Corden, J Sinha, R Jones. Dry cotton or flocked respiratory swabs as a simple collection technique for the molecular detection of respiratory viruses using real-time NASBA. *J of Virol Methods* ۱۵۳ (۲۰۰۸) ۸۴-۸۹.
۲. Boixeda R, Rabella N, Sauca G, Delgado M, Martínez-Costa X, Mauri M, Vicente V, Palomera E, Serra-Prat M, Capdevila JA. Microbiological study of patients hospitalized for acute exacerbation of chronic obstructive pulmonary disease (AE-COPD) and the usefulness of analytical and clinical parameters in its identification (VIRAE study). *Int J Chron Obstruct Pulmon Dis*. ۲۰۱۲; ۷:۳۲۷-۳۵.
۳. Debyle C^۱, Bulkow L, Miernyk K, Chikoyak L, Hummel KB, Hennessy T, Singleton. Comparison of nasopharyngeal flocked swabs and nasopharyngeal wash collection methods for respiratory virus detection in hospitalized children using real-time polymerase chain reaction. *J Virol Methods*. ۲۰۱۲ Oct; ۱۸۵(۱):۸۹-۹۳.
۴. Falsey AR, Formica MA, Walsh EE. Yield of sputum for viral detection by reverse transcriptase PCR in adults hospitalized with respiratory illness. *J Clin Microbiol*. ۲۰۱۲ Jan; ۵۰(۱):۲۱-۴.
۵. Bartlett JG. Diagnostic tests for agents of community-acquired pneumonia. *Clin Infect Dis*. ۲۰۱۱ May; ۵۲ Suppl ۴:S۲۹۶-۳۰۴.
۶. Pongthanapisith V, Sukasem C, Premchaiporn K, Srichantaratsamee C, Chantratita W. Clinical performance of three rapid diagnostic tests for influenza virus in nasopharyngeal specimens to detect novel swine-origin influenza viruses. *Infection*. ۲۰۱۱ Apr; ۳۹(۲):۱۰۵-۱۱.
۷. Piralla A, Giardina F, Rovida F, Campanini G, Baldanti F. Cellular DNA quantification in respiratory samples for the normalization of viral load: a real need? *J Clin Virol*. ۲۰۱۸ Oct; ۱۰۷:۶-۱۰.

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